



## E.coli and total coliforms, pk/25

Nr. produs:	2495525
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Presence/Absence Broth with MUG ampules provide an easy-to-use test, ideal for screening drinking water samples for the presence of Total Coliforms and E. coli bacteria. Simply combine the P/A broth with your 100 mL sample and incubate for 24 to 48 hours at 35 °C. A colour change to yellow or yellow brown indicates presence of coliform bacteria.

**Complete Solutions for Testing Recreational Water**

**Complete Solutions for Testing Wastewater**

**Complete Solutions for Testing Drinking Water**

---

### Specificații

Cantitate:	25 /pk
Material:	Glass