



## Ampules, m-ColiBlue24 broth, glass , pk/20

**Nr. produs:** 2608420

**RON Preț (fără TVA):** Contact

**Expediere în 2 săptămâni**

m-ColiBlue24 Broth Ampules; Glass

m-ColiBlue24 Broth Ampules; Glass

In 24 hours, Hach's m-ColiBlue24® Broth give you simultaneous results for total coliforms and E. coli. Hach's ready-to-use media eliminates measuring, mixing, and autoclaving steps necessary to prepare dehydrated media. Glass ampules (26084-20), plastic ampules (260845-00), 100mL bottles (26084-42) of m-ColiBlue24® and agar plates (28052-15), afford maximum shelf life and ease of use for pre-measured and prepared media. Enumerate total coliforms and E. coli on one petri dish. Differentiate colonies easily - red and blue indicate total coliforms and blue specifies E. coli. Get optimal recovery of stressed and injured organisms.

Use m-ColiBlue24® Broth to accurately monitor and evaluate total coliform and E. coli in drinking water; wastewater; bottled water; beverages; surface, ground, and well water; and ultrapure, chemical processing and pharmaceutical processing waters.

m-ColiBlue24 should only be used with pads manufactured by the following companies/brands: Pall Corporation (e.g., Sterile Petri Dishes with Absorbant Pads), Private Label/OEM products manufactured by Pall, Gelman, Sartorius (cellulose and glass fiber pads), Sartorius BioSart products with the integrated Petri dish. No other pad is compatible with the product, and the use of any pad other than those recommended may result in unreliable, inaccurate data.

With ready-to-use m-ColiBlue24® Broth, you can:

**Read and confirm results in 24 hours**

**Get superior sensitivity - 1 CFU/100mL**

**Conduct emergency testing and routine monitoring**

**Eliminate special equipment and ultraviolet lamp**

**Minimize background growth of non-coliform bacteria**

---

### Specificații

Cantitate:	20 /pk
Material:	Glass
Media Type:	m-ColiBlue24 Broth