



Reagent set, aluminum, 50 mL sample, ECR method

Nr. produs:

2443400

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



For the determination of Aluminium, Eriochrome Cyanine R Method. For 25ml sample.

For the determination of Aluminium, Eriochrome Cyanine R Method. For 25ml sample.

For determination of Aluminum by the Eriochrome Cyanine R (ECR) Method. Hach Method 8326. Sample Size: 50 mL Range: 0.002 - 0.220 mg/L. Reagent set contains ECR Reagent Powder Pillows, Hexamethylenetetramine Buffer Reagent Powder Pillows and ECR Masking Reagent Solution for approximately 100 tests.

Specificații

Domeniu de măsurare:	0 - 0,220 mg/L Al
Instrument:	DR/3000, DR/2000
Numar metoda:	8326
Numar teste:	100
Nume metoda:	Eriochrome Cyanine R