



Mercuric nitrate cartridge for digital titrator, 0.2256 N, 1 piece

Nr. produs:

1439301

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Ready-to-use titrant solution for the analysis of Chloride (Mercuric Nitrate Method)

Ready-to-use titrant solution for the analysis of Chloride (Mercuric Nitrate Method)

For chloride determination by Digital Titrator titration with mercuric nitrate, using diphenylcarbazone indicator. Hach method 8206. Digits = mg/L Cl⁻ in 100 mL sample. Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	10 - 160 mg/L Cl ⁻
Numar metoda:	8206
Numar teste:	40 - 100
Nume metoda:	Mercuric Nitrate
Parametru:	Chloride
Platforma:	Titrator digital
Titrant:	0.2256 N Hg(NO ₃) ₂