



Reagent set, lead, dithizone method, (no air shipments)

Nr. produs:

2243100

Produse periculoase

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.



For the determination of Lead, Dithizone Method.

For the determination of Lead, Dithizone Method.

For determination of lead by the Dithizone extraction method. USEPA accepted for reporting for wastewater analysis using Hach Method 8033. Range: 3 to 300 µg/L Pb. Set Includes reagents for 100 tests. For air freight shipments, order 2243101

Specificații

Domeniu de măsurare:	3 - 300 µg/L Pb
EPA compliant:	Da
Instrument:	All except DR/800, PCII
Numar metoda:	8033
Numar teste:	60 - 100
Nume metoda:	Dithizone
Shipping Method:	Lead Reagent set: No Air Freight