



Be Right™



Reagent set, dissolved oxygen, buret, 300 mL sample

Nr. produs: 2351400
RON Preț (fără TVA): Contact
Apelare pentru data de expediere

Produse periculoase



Reagents for determination of Dissolved Oxygen by Winkler method titration

Reagents for determination of Dissolved Oxygen by Winkler method titration

For determination of Dissolved Oxygen by the Winkler Method (Azide Modification) and buret titration. Range: up to 10 mg/L as DO. Sample Size: 300-mL. Contains reagents for 100 tests.

Specificații

Domeniu de măsurare:	1 - > 10 mg/L O ₂
EPA compliant:	Da
Numar metoda:	8229
Numar teste:	100