



Glass fiber filter, pore size 1.5 µm, 47 mm diameter, 100 pieces

Nr. produs:	253000
RON Preț (fără TVA):	Contact
Disponibil	

for TSS (Total Suspended Solids)

for TSS (Total Suspended Solids)

TSS Filters are recommended for use in Standard Methods 2540 C and 2540 D. Optimized surface characteristics of this binderless borosilicate glass microfiber filter eliminates pin-holing and fiber-shedding to ensure consistent, accurate results. Fast flow rate, fine porosity and 1.5 µm particle retention makes the TSS Filter the preferred choice of environmental chemists for total suspended solids analysis of water and wastewater.

Specificații

Cantitate:	100 /pk
Diametru:	47 mm
Dimensiune por:	1.5 µm
Material:	Glass Microfiber