



## Bag, Whirl-Pak\* 5-1/2X9, pk/10

Nr. produs: 1437297

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Complies with Standard Methods and USEPA for drinking water sample collection. Made from poly resins approved by USDA and FDA for food contact. Sterile (certified), disposable, transparent polyethylene bags with trusted "Whirl-Pak" closure. Suitable for liquid or solid sample collection.

---

### Specificații

Cantitate:	10
Capacitate:	720 mL