



Hardness, iron, and pH test kit, model HA-62A

Nr. produs:

183701

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



Model HA-62A

Model HA-62A

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specificații

Aplicatie:	Water Quality
Domeniu de măsurare:	6.6 - 8.4 pH units
Model:	HA-62A
Numar teste:	100
Parametru:	Hardness, Iron, and pH (6.6-8.4)
Platforma:	Disc culoare