



Test kit, chlorine, hardness, iron & pH , model CN-39WR

Nr. produs: 223002
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Produse periculoase



Chlorine, Hardness, Iron, and pH Test Kit combines the most popular water treatment test into one compact kit. The chlorine test is USEPA-accepted for drinking water testing if approved by individual states.

Includes frequently measured water quality parameters in one rugged test kit

Uses simple, economical drop count titration and color disc methods

Tests chlorine (0-3.4 mg/L), hardness (1-20 gpg), iron (0-4 mg/L), and pH (4-10).

Specificații

Aplicatie:	Water Quality
Domeniu de măsurare:	0 - 3.4 mg/L
Model:	CN-39WR
Numar teste:	100
Parametru:	Chlorine, Hardness, Iron, and pH
Platforma:	Disc culoare