



Be Right™



Bottle, Wash, Nalgene, Narrow Mouth, 500 mL, Red Cap/ Stem, 6/pk

Nr. produs: 2928204

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Narrow mouth polypropylene screw closure and stem feature one-piece construction to ensure a leak proof seal. Tip can be cut back to increase flow. 500 mL (16.8 oz. Capacity, 28 mm mouth). Choose from 4 cap/stem colors.

Specificații

Cantitate:	6 pk
Capacitate:	500 mL
Culoare:	Incolor
Material:	LDPE
Stem:	Red Cap/Stem
Unitate:	6 pk