



Beaker, PP, low form, 100 mL

Nr. produs: 108042

RON Preț (fără TVA): Contact

Disponibil

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specificații

Cantitate:	1 each
Capacitate:	100 mL
Diametru exterior :	63.0 mm with spout
Material:	Polypropylene, Low Form
Material corp:	Polypropylene
Unitate:	1 each