



Beaker, polypropylene, low form, 50 mL

Nr. produs:	108041
RON Preț (fără TVA):	Contact
Disponibil	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specificații

Cantitate:	1 each
Capacitate:	50 mL
Diametru exterior :	51.1 mm with spout
Material:	Polypropylene, Low Form
Material corp:	Polypropylene
Unitate:	1 each