



## Test kit, acidity, model AC-DT

**Nr. produs:** 2064000  
**RON Preț (fără TVA):** Contact  
 Nu este afișată data expedierii

Produse periculoase



Acidity is usually caused by the presence of mineral acids, salts of strong acids, or free carbon dioxide. In surface or ground waters, acidity may form natural substances or industrial pollution. Excessive levels cause corrosion and can be detrimental to fish. Utilizing the Hach Digital Titrator, the Model AC-DT kit is used to determine acidity concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-4,000 mg/L. Contains approximately 100 tests.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

### Specificații

Cei mai mici pași de incrementare:	0.1-1.10
Domeniu de măsurare:	10 - 4,000 mg/L CaCO <sub>3</sub>
Model test kit:	AC-DT
Nota de subsol:	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. 1EPA-accepted for wastewater. Check with your regulatory agency.
Numar teste:	100
Nume metoda:	Digital Titrator/Sodium hydroxide
Parametru:	Acidity - (P or Total) as CaCO <sub>3</sub>
Platforma:	Titration digital
Tipul carcasei:	F