



## Bottle Carrier, Linear

**Nr. produs:** 2094300

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Linear, low-density polyethylene. Securely holds up to six standard 0.5 L glass bottles, pint-size reagent bottles, and most other 500 mL glass or plastic bottles—making breakage almost impossible. Lightweight, chemical-resistant, and easy to clean. Dimensions: 34.9 x 24.8 x 27.6 cm (14 x 10 x 11"). Maximum bottle size it will hold: 73 mm in diameter, 178 mm in height.

---

### Specificații

|            |           |
|------------|-----------|
| Adâncime:  | 11 Inches |
| Cantitate: | 1 each    |
| Latime:    | 10 Inches |
| Lungime:   | 349 mm    |