



Agar R2A, dehydrated 500 g

Nr. produs: 2281034

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Low nutrient medium for enumeration and cultivation of bacteria from potable water. Used for pour or spread plate and MF methods; described in the 16th edition of Standard Methods. Recommended for predicting nitrification more accurately in potable systems using chloramines as a secondary disinfectant.

Specificații

| | |
|-------------|----------|
| Cantitate: | 500 g |
| Media Type: | R2A Agar |