



Test kit, surface water

Nr. produs:

2559800

RON Preț (fără TVA):

Contact

Nu este afișată data expirării

Produse periculoase



The ideal tool for source water monitoring and wastewater treatment plants, this kit uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, nitrate, dissolved oxygen, phosphorus, and temperature. The rugged, lightweight kit includes the convenient pH Pocket Pal Tester and all necessary apparatus and reagents for testing.

Ideal tool for both industrial and municipal water treatment plants

Uses drop count titration and color disc methods to provide economical, easy-to-use tests for ammonia, chlorine, pH, titrate, dissolved oxygen, phosphorus, and temperature

Rugged, lightweight kit includes the convenient pH Pocket Pal[™] Tester and all necessary apparatus and reagents

Specificații

Aplicatie:	Surface Water
Numar teste:	100
Platforma:	Disc culoare