



## Cylinder, graduated, mixing, certified, 100 mL

Nr. produs: 2636342

RON Preț (fără TVA): Contact

Nu este afișată data expirării

This certified, Class A glass mixing cylinder is calibrated “to contain” according to ASTM E542 & E1272. Serialized and traceable to NIST standards (includes Certificate of Identification and Capacity). Single metric scale in durable white markings with funnel-top. Glass penny head stopper included.

---

### Specificații

Cantitate:	1 each
Capacitate:	100 mL
Certificări:	Class A
Material corp:	Glass
Metoda de măsurare:	To Contain
Tip:	Simplu Metric, Mixing