



Buret, Class A, Glass, 25 mL

Nr. produs: 2636540

RON Preț (fără TVA): Contact

Nu este afișată data expedierii

Glass, funnel-fill buret. Certified Class A (ASTM Specification E287) with permanent serial number and certificate for NIST traceability. PTFE stopcock does not require grease. Dust cap included. 10 mL buret meets requirements for potentiometric titrations (ASTM D664).

Specificații

Capacitate:	25 mL
Caracteristici speciale:	Buret
Material corp:	Class A Glass