



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Nr. produs: 2936500

RON Preț (fără TVA): Contact

Nu este afișată data expirării

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 47 mm

Specificații

Cantitate:	600
Diametru:	47 mm mm
Dimensiune por:	0.45 µm
Material:	Mixed-Cellulose Ester