



Fixed Volume, Micropipette, 10 µL

Nr. produs: 2946002

RON Preț (fără TVA): Contact

Nu este afișată data expirării

Micropipette features a two-stage plunger dispenser, which ensures high accuracy and precision in pipetting techniques and facilitates last drop dispensing. Unique tip cone accepts both micro and ultra micro pipet tips. Each micropipette is fully autoclavable and colour-coded for easy identification. Includes one pipet tip. Additional pipet tips sold separately.

Specificații

Acuratete:	1,0 %
Cantitate:	1 each
Capacitate:	10 µL