



Beaker, Griffin, Low Form, Glass, 600 mL

Nr. produs: 50052
RON Preț (fără TVA): Contact
Nu este afișată data expirării

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specificații

Capacitate:	600 mL
Material:	Glass, Griffin, Low Form
Material corp:	Glass
Unitate:	1 each