



Pachet de pornire TOC/TNb QP1680

Nr. produs: **SMKIT500000**

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Specificații

Conținutul cutiei:

- 1x Cotton swabs, Cotton/Plastic, 25 pcs
- 1x Furnace injection port, PTFE
- 1x Filter, Plastic, 5 µm and disposable, 10 pcs
- 1x Stirrer bar, PTFE, 15 pcs
- 5x O-ring Injection and IC port
- 1x Allen key, 4 mm
- 1x Leak detect assembly
- 1x Pressure test tool (up to 1 bar)
- 1x Septa Aluminium, 1000 pcs
- 1x Tweezer, 250 mm
- 1x Needle, Stainless steel
- 1x Use tool, O-ring replacement
- 1x Service kit needle guide

Conținutul cutiei

1x Cotton swabs, Cotton/Plastic, 25 pcs
1x Furnace injection port, PTFE
1x Filter, Plastic, 5 µm and disposable, 10 pcs
1x Stirrer bar, PTFE, 15 pcs
5x O-ring Injection and IC port
1x Allen key, 4 mm
1x Leak detect assembly
1x Pressure test tool (up to 1 bar)
1x Septa Aluminium, 1000 pcs
1x Tweezer, 250 mm
1x Needle, Stainless steel
1x Use tool, O-ring replacement
1x Service kit needle guide