



## Acid solution for EZ4004 Total Alkalinity Analyser, 0.2 N HCl (100-5000 mg/L Range), 2 L

Nr. produs:

APPC4004-04-EU **Produse periculoase**

RON Preț (fără TVA):

Contact

Nu este afișată data expirării



### Ready-to-use Reagents for EZ Series Online Analysers

#### Ready-to-use

Hach EZ reagents are ready-to use with EZ Series Online Analysers. Since there is no additional preparation needed this eliminates an additional error source when using the EZ instruments and reagents system.

#### Confidence in your results

Reduce the sources of errors in your workflow to increase the confidence in your results.

#### Less complexity

Using ready-to-use EZ reagent kits from Hach can reduce the number of handling steps around your analytical equipment significantly which allows less hands-on-time and complexity.

---

### Specificații

Cantitate:	2 L
Concentrație:	0.2 N HCl
Parametru:	Total Alkalinity