



Sedimenter for AISE, NISE and AN-ISE sc

Nr. produs:	LZY771
RON Preț (fără TVA):	Contact
Expediere în 1 săptămână	

The sedimenter can be used as a bypass solution. The sample can be let through the sedimenter without prefiltering.

Specificații

Greutate:	Sedimenter incl. wall plate only: 14.7 kg
Instrument:	Suitable for all AISE, NISE and AN-ISE sensors
Limbi manual:	Non verbal instruction sheet
Lungime:	Intallation plate for wall mount: 800 mm x 495 mm
Material:	PVC
Modalitate de montaj:	Wall installation (indoor only).