



Be Right™



## Sistem montare Filtrax/Sonda filtrare, brida 10 cm, bara INOX 2m

**Nr. produs:** LZY714.99.42020  
**RON Preț (fără TVA):** Contact  
**Disponibil**

---

### Specificații

Diametru:	40 mm
Greutate:	Lifting part only: 5.5 kg
Instrument:	suitable for:  FILTERPROBE sc;  FILTRAX;  SIGMATAX .
Limbi manual:	non verbal instruction sheet
Lungime:	Pole: 2000 mm
Material:	Wetted parts:  pole SS 1.4301;  seals EPDM;  screw set SS.