



Ammonia Ionic Strength Adjustor (ISA) Powder Pillows, pk/100

Nr. produs:

4447169

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Ammonia Ionic Strength Adjustor Buffer for use when making ammonia determinations by electrode methods. Ideal for determining low ammonia levels.

Ammonia Ionic Strength Adjustor Buffer for use when making ammonia determinations by electrode methods. Ideal for determining low ammonia levels.

For use in ammonia determinations by electrode methods. Powder form does not dilute samples. Ideal for determining low ammonia levels. Pack of 100 Powder Pillows

Specificații

Cantitate:	100 /pk
Concentratie:	ISA
Package settings:	Powder Pillows
Parametru:	Ammonium