



Free and total chlorine test strips, 0 - 10 mg/L, 50 tests

Nr. produs: 2745050
RON Preț (fără TVA): Contact
Expediere în 2 săptămâni

In the field, technicians can test many samples in only a few minutes, and make immediate evaluations on-site. No measuring, set-up, clean-up, or chemical handling are necessary. Lab technicians also use the test strips for pre-test screening-to detect the presence of materials that might interfere with testing. Steps: (mg/L) 0, 0.5, 1.0, 2.0, 4.0, 10.0.

Specificații

Cei mai mici pasi de incrementare:	Steps: 0, 0.5, 1, 2, 4, 10
Domeniu de măsurare:	0 - 10 mg/L Cl_2
Model test kit:	TS
Nota de subsol:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Parametru:	Free and Total Chlorine
Platforma:	Test Strip
Unitate:	50 tests