



Be Right™



pH Test strips, 4 - 9 pH units, 50 tests

Nr. produs: 2745650
RON Preț (fără TVA): Contact
Disponibil

Most natural waters range from pH 4 to pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range. Contains 50 tests.

Easy to use, disposable, and inexpensive

Great for quick evaluations of pH levels in the lab or in the field

Specificații

| | |
|------------------------------------|----------------------------|
| Cei mai mici pasi de incrementare: | Steps: 4, 5, 6, 7, 8, 9 pH |
| Domeniu de măsurare: | 4 - 9 pH units |
| Model test kit: | Test Strips |
| Parametru: | pH, mid range |
| Platforma: | Test Strip |
| Unitate: | 50 tests |