



**Be Right™**



## Test strips pH , 4-9 pH units

**Nr. produs:** 2745650

RON Preț (fără TVA): Contact

Disponibil

Most natural waters range from pH 4 to pH 9, but commonly are above pH 7 because of carbonates and bicarbonates (alkalinity). The pH of water used in industry, boiler feed water, and swimming pools usually is kept within a narrow range. Contains 50 tests.

**Easy to use, disposable, and inexpensive**

**Great for quick evaluations of pH levels in the lab or in the field**

---

### Specificații

Cei mai mici pasi de incrementare:	Steps: 4, 5, 6, 7, 8, 9 pH
Domeniu de măsurare:	4 - 9 pH units
Model test kit:	Test Strips
Parametru:	pH, mid range
Platforma:	Test Strip
Unitate:	50 tests