



Be Right™



Filter, glass fiber, diameter 55 mm, preweighed standard methods TSS procedure, 100/pk

Nr. produs: 2546120

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Individually preweighed filter papers, preweighed in compliance with Standard Methods Total Suspended Solids procedure. Each in a planchet with tare weight listed. Each tare weight meets Standard Methods specifications.

Specificații

Cantitate:	100 /pk
Diametru:	55 mm
Dimensiune por:	1.5 μ m
Material:	Glass Microfiber