



Soxhlet thimbles, cellulose, 33 x 80 mm, 25 pk

Nr. produs:

2966700

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

Xenosep® extraction thimbles are ideal for determining the concentration of nonvolatile and semivolatile organic compounds from environmental solid waste matrices, sludge, soil, air, and water.

The thimbles are manufactured from high purity alpha cellulose cotton linters which provide excellent mechanical strength, stability, and retention with very low/no background extractables in a disposable, sample holder format.

Specificații

Cantitate:	25 /pk
Dimensiune por:	10 µm
Material:	Cellulose