



Be Right™



Dispenser, variable volume, digital, 1.0-10.0 mL

Nr. produs: 2563140

RON Preț (fără TVA): Contact

Apelare pentru data de expediere

Variable volume bottle-top dispensers screw directly onto reagent bottles to deliver solution volumes within tightest tolerances. Fixed selector steps permanently calibrated at factory.

Variable volume bottle-top dispensers screw directly onto reagent bottles to deliver solution volumes within tightest tolerances. Fixed selector steps permanently calibrated at factory.

Attach dispenser directly to reagent bottle for adjustable dispensation with a single stroke. Inert hardware designed to handle alkaline, acidic and organic solutions. Autoclavable (121°C). Dispensers fits bottles with a 33 mm screw-neck finish (supplied with adapters to fit 28, 38 and 45 mm finishes). Dispenser rotates 360°.

Specificații

Cantitate: 1 each