



## Filter Holder, Magnetic

Nr. produs:	1352900
RON Preț (fără TVA):	Contact
Expédiere în 2 săptămâni	

Clear polysulfone funnel with imprinted graduations. Magnetic coupling of funnel to stem for fast filter changes. Used with 47 mm membrane filter for coliform testing and prefiltration of many fluids.

Clear polysulfone funnel with imprinted graduations. Magnetic coupling of funnel to stem for fast filter changes. Used with 47 mm membrane filter for coliform testing and prefiltration of many fluids.

Fits standard manifold systems with 8 rubber stopper connections. Capacity 300 ml (10 oz). Measures 23 x 7 cm (9 x 2.75"). Autoclavable.

---

### Specificații

Cantitate:	1
------------	---