



## Tube, culture, w/ screw cap, 25 x 150 mm, 10/pk

Nr. produs: 2096250

RON Preț (fără TVA): Contact

Nu este afișată data expirării

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

The cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

Cap has a firmly bonded rubber liner, providing a sure and highly inert sealing surface.

---

### Specificații

Cantitate:	10 /pk
Capacitate:	Capacity to Neck 50 mL
Caracteristici speciale:	Tube, 25 x 150 mm
Diametru:	25 mm
Dimensiuni (A x L):	25 mm x 150 mm
Lungime:	150 mm