



Be Right™



Bottle, sampling, presterilized, round, 50/pk

Nr. produs: 2495050

RON Preț (fără TVA): Contact

Expediere în 2 săptămâni

Hach microbiological chemicals

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen(TM) Medium.

Specificații

Cantitate: 50/pk