



## Paddle tester, total aerobic bacteria, yeast and mold, 10 pieces

<b>Nr. produs:</b>	<b>2610810</b>
RON Preț (fără TVA):	Contact
Expediere în 2 săptămâni	

Hach Paddle Testers allow fast, easy screening of water samples, solid surfaces, and nonviscous liquids where unsanitary conditions may exist. The paddle is a double-sided slide attached to the vial cap. Each side of the slide is used to perform a separate test. The slide has a molded-in grid for easy colony counting without opening the optically clear, leakproof vial. Paddle Tester, Total Aerobic Bacteria/Yeast and Mold (Tryptone Glucose Extract Agar/Sabouraud Dextrose Agar), pk/10.

### Save time and simplify lab and field testing with Hach prepared P/A media

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Report Total Coliform values with EPA accepted Presence-Absence bottles and media

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Simplify routine testing with self-contained, easy-to-use Paddle Testers and BART products

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Specificații

Cantitate:	10 /pk
Media Type:	Total Aerobic Bacteria, Yeast and Mold