



C501X Cable SubD-9f/2m/SubD-9f

Nr. produs:

A95X501

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

RS232 Serial Connection Cable, SubD-9

RS232 Serial Connection Cable, SubD-9

The CCTO is a pyrex glass body cell designed for testing corrosion on large flat samples. The shape of the cell combined to the large area of the working electrode and the stainless steel counter electrode make it particularly suitable for testing low conductivity media and/or highly resistive samples. The setup includes a calomel reference electrode mounted on a vertically adjustable Luggin capillary designed to minimize ohmic drop.