



Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe, FerroVer

Nr. produs:

85428

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Powder reagent in sealed package, for 25mL sample.

Powder reagent in sealed package, for 25mL sample.

FerroVer Iron Reagent: For the determination of total soluble Iron including all ferric (Fe³⁺) and ferrous (Fe²⁺) Iron as well as complexed Iron compounds (Fe-EDTA). Perfect for many applications and a wide range of concentrations. Determines Iron in the presence of Copper.

Specificații

Cantitate:	1000 /pk
Domeniu de măsurare:	0,02 - 3,00 mg/L Fe
EPA compliant:	Da
Numar metoda:	8008
Numar teste:	1000
Nume metoda:	FerroVer
Platforma:	Pernu#e cu pudră