



FerroVer Iron Reagent Powder Pillows, 0.02-3.00mg/L Fe

Nr. produs:

85499

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Powder reagent in sealed packages for the determination of Iron by the 1,10 phenantroline method. pk/100. For 25 mL samples.

Specificații

Cantitate:	100 /pk
Condiții de depozitare:	10 °C - 25 °C
Descriere:	FerroVer Iron Reagent Powder Pillows, 25 ml
Domeniu de măsurare:	0,02 - 3,00 mg/L Fe
EPA compliant:	Da
Instrument:	DR/3000, DR/2000
Necesitate digestie:	For total iron
Numar metoda:	8008
Numar teste:	100
Nume metoda:	FerroVer
Parametru:	Iron
Platforma:	Pernu#e cu pudră