



ZincoVer 5 Zinc Reagent Powder Pillows

Nr. produs:

2106669

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Powder reagents in sealed package for the determination of Zinc, Zincon Method. For 20mL sample.

Powder reagents in sealed package for the determination of Zinc, Zincon Method. For 20mL sample.

Specificații

Cantitate:	100/pk
Domeniu de măsurare:	0,01 - 3,00 Zn
EPA compliant:	Da
Numar metoda:	8009
Numar teste:	100
Nume metoda:	Zincon
Parametru:	Zinc
Platforma:	Pernu#e cu pudră