



CuVer 2 Copper Reagent Powder Pillows

Nr. produs:

2188299

Produse periculoase

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni



For the determination of copper, 0.04-5.00mg/L Cu, in the presence of high levels of Aluminum, Calcium, and Iron. For 25 mL sample. Pack of 100.

For the determination of copper, 0.04-5.00mg/L Cu, in the presence of high levels of Aluminum, Calcium, and Iron. For 25 mL sample. Pack of 100.

Specificații

Cantitate:	100 /pk
Numar teste:	100
Nume metoda:	Bichinchoninate
Parametru:	Copper
Platforma:	Pernu#e cu pudră