



Digital titrator cartridge, sulfuric acid, 1.600 N, 1 piece

Nr. produs:

1438901

RON Preț (fără TVA):

Contact

Disponibil

Produse periculoase



Ready-to-use titrant solution for acid-base and alkalinity titrations.

Specificații

Domeniu de măsurare:	100 - 4000 mg/L CaCO_3
Domeniu de masurare :	1 - 10 meq/L Base
EPA compliant:	Da
Method Number 1 Chemistry:	8203
Numar teste:	40 - 100
Nume metoda:	Verde de bromcresol/rojo de metilo
Parametru:	Alkalinity, Total
Platforma:	Titratore digital
Titrant:	1.600 N H_2SO_4