



TitraVer EDTA Digital Titrator Cartridge, 0.800 M

Nr. produs:

1439901

RON Preț (fără TVA):
Disponibil

Contact

Produse periculoase



Ready-to-use titrant solution for the determination of Calcium and Total Hardness, 100-4000mg/L CaCO₃.

Specificații

Domeniu de măsurare:	100 - 4000 mg/L CaCO ₃
Numar metoda:	8213
Numar teste:	40 - 100
Nume metoda:	EDTA (ManVer® 2)
Parametru:	Dureza, total, secuencial
Platforma:	Titrator digital
Titrant:	0.800 M EDTA