



## EDTA Digital Titrator Cartridge, 0.714 M

**Nr. produs:**

**1495901**

**Produse periculoase**

RON Preț (fără TVA):

Contact

Disponibil



Ready-to-use titrant solution for the analysis of Total Hardness, 10-100+ G.d.H. CaCO<sub>3</sub>.

### Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

### Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

### Versatile

More than 50 parameter/range combinations are available

### Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

---

## Specificații

|                      |                   |
|----------------------|-------------------|
| Domeniu de măsurare: | 10 - 200 G.d.h.   |
| Numar metoda:        | 8204              |
| Numar teste:         | 40 - 100          |
| Nume metoda:         | EDTA (CalVer2)    |
| Parametru:           | Hardness, Calcium |
| Platforma:           | Titratore digital |
| Titrant:             | 0.714 M EDTA      |