



SulfaVer 4 Sulfate Reagent Powder Pillows, 2-70mg/L SO₄

Nr. produs:

1206599

RON Preț (fără TVA):

Contact

Disponibil

Produse periculoase



Powder reagents in sealed package for determination of sulfate. For 25 mL samples. pk/100

Specificații

Cantitate:	100 /pk
Descriere:	SulfaVer 4 Sulfate Reagent Powder Pillow, 25 ml
Domeniu de măsurare:	2 - 70 mg/L SO ₄
EPA compliant:	Da
Instrument:	DR/4000, DR/3000, DR/2010, DR/2000
Nota de subsol:	*Generally, the DR 5000 performs the same methods as the DR/4000. The DR 2800 performs the same methods as the DR 2700, DR/2500, and DR/2400. Some exceptions do exist.
Numar metoda:	8051
Numar teste:	100
Nume metoda:	SulfaVer 4
Parametru:	Sulfate
Platforma:	Pernu#e cu pudră
Reference method:	USEPA 375.4