



Iodate-Iodide Digital Titrator Cartridge, 0.3998 N

Nr. produs:

1496101

RON Preț (fără TVA):

Contact

Expediere în 2 săptămâni

Produse periculoase



Ready-to-use titrant solution for the analysis of Sulfite, 4->400mg/L SO₃.

Ready-to-use titrant solution for the analysis of Sulfite, 4->400mg/L SO₃.

Iodate-Iodide Standard Solution. For determination of sulfite by Digital Titrator iodometric titration. Digits = mg/L SO₃ in 20 mL sample.
Digits x 0.0005 = meq. Approximately 40-100 tests.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specificații

Domeniu de măsurare:	0 - 800 mg/L SO ₃ ²⁻
Nume metoda:	Iodate-Iodide
Parametru:	Sulfite
Platforma:	Titrator digital
Titrant:	0.3998 N KIO ₃ -KI