



Sension+ DO6 Portable Dissolved Oxygen Meter Field Kit with Sensor

Nr. produs:

LPV4552.98.0002

ARTICOL ÎNVECHIT

Acest articol nu mai este disponibil.

Complete mobile measurement stations.

Kit contains portable DO meter plus polarographic DO probe 5130, field case and accessories. For general aqueous applications.

Kit contains portable DO meter plus polarographic DO probe 5130, field case and accessories. For general aqueous applications.

The basic Dissolved Oxygen field kit is an all-in-one portable laboratory that makes general water quality meter testing in the field fast and simple. The field kit is designed to be used in a wide variety of applications and comes complete with everything needed to start testing.

Simple and fast measurements

Ergonomic design for real one hand operation

IP67 protected

Precise and cost-efficient system.

Specificații

Acuratete:	-
Acuratete temperatura:	0,2 °C (± 1 cifră)
Acuratetea conductivitatii:	-
Autorecunoasterea solutiei buffer:	-
Calibrare electrod oxigen diolvat:	1 to 2 point, saturated air, 0% DO calibration chamber delivered with kit
Capacitate de stocare rezultate interne:	-
Caracteristici GLP:	-
Cerință acumulator:	3 AA
Certificarea conformității:	CE.WEEE
Conditii de mediu: temperatura:	0 - 50 °C (32 - 122 °F)
Condiții de mediu: umiditate relativă:	< 80 % relative humidity (non-condensing)
Dimensiuni (Înălțime x Lățime x Adâncime):	38 mm x 74 mm x 186 mm
Electrod inclus:	5130 Electrode
Evaluare IP carcasa:	IP67
Garanție:	36 luni

Greutate:	0,3 kg
Ieşire:	-
Iluminare de fundal:	Yes
Inaltime (global):	-
Interfață operare:	Keypad
Interfata: limbi:	Icon based, no language needed
Intrari:	MP-5
Kit?:	Da
Limbi interfață utilizator:	-
Limbi manual:	Multi-lingual
Masa de lucru:	No
Masurare oxigen dizolvat:	0.00 - 19.99 mg/L and 20.0 - 22.0 mg/L (25°C)
Masurare salinitate:	-
Masurare TDS:	-
Masurare temperatura:	-20.0 - 150.0 °C (-4.0 - 302.0 °F)
Masurarea conductivitatii:	-
Masurarea conductivitatii la citire stabila:	-
Masurarea conductivitatii: corectia temperaturii:	-
Masurarea presiunii barometrice:	Manual introduction
Measuring range:	-
Metoda de măsurare:	Auto stabilisation, continuous
Model:	DO6
Monitor:	Backlit LCD with pictograms
Operare AC si USB:	No
Parametru:	Dissolved Oxygen, Temperature
PC data transfer software :	-
pH Accuracy:	-
pH Electrode calibration:	-
pH Measurement Range:	-
Pressure compensation:	Manual input
Rezistenta la apa:	IP67
Rezolutia conductivitatii:	-
Rezoluție:	0.01 mg/L
Rezolutie mV:	-
Rezolutie oxigen dizolvat:	0.01 mg/L depending on measurement range
Rezolutie pH:	-
Rezolutie TDS:	-
Rezolutie temperatura:	0,1 °C
Salinity resolution:	-
Simultaneous measurements:	Dissolved Oxygen, Temperature
Temperature compensation:	Automatic or manual
Tip display:	128 x 64 pixel, simultaneous readings

Accesorii necesare

- Soluție de umplere, internă, sondă polarografică pentru oxigen dizolvat, 30 mL (Item 2759123)
- Soluție de umplere, internă, sondă polarografică pentru oxigen dizolvat, 50 mL (Item 2759126)