



CDC866T General-purpose conductivity cell, 4-pole, Temp.

Nr. produs: E61M015
RON Preț (fără TVA): Contact
Nu este afișată data expedierii

4-pole conductivity electrode with MAB6 plug and epoxy body. Cell constant 1cm^{-1} . Platinised, with temperature sensor. Max. temperature 80°C , min. immersion depth 35mm.

Specificații

Caracteristici speciale:	Conductivity electrode with MAB6 plug and epoxy body.
Conductivity Cell:	CDC866T
Electrod:	4-pole, platinised
Garanție:	12 luni