



Glass sample cell for lab turbidimeters, 30 mL, 6 pieces with caps

Nr. produs: 2084900

RON Preț (fără TVA): Contact
Disponibil

For use with 2100N, 2100AN, and Lab RATIO turbidimeters with caps

Specificații

Greutate:	0,2 kg
Material:	Glass
Nota de subsol:	For use with Hach Ratio 18900, 2100N and 2100AN Laboratory Turbidimeters.
Scop alimentare:	6 vials with caps
Volum probă:	30 mL